

## Product datasheet for **MR222965L4V**

### Glp1r (NM\_021332) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Glp1r (NM_021332) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Glp1r
Synonyms:	GLP-1R; GLP1Rc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021332
ORF Size:	1392 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222965).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021332.2</a> , <a href="#">NP_067307.2</a>
RefSeq Size:	1480 bp
RefSeq ORF:	1392 bp
Locus ID:	14652
UniProt ID:	<a href="#">O35659</a>
Cytogenetics:	17 15.8 cM



[View online »](#)

**Gene Summary:**

G-protein coupled receptor for glucagon-like peptide 1 (GLP-1) (PubMed:9568699). Ligand binding triggers activation of a signaling cascade that leads to the activation of adenylyl cyclase and increased intracellular cAMP levels (By similarity). Plays a role in regulating insulin secretion in response to GLP-1 (PubMed:9568699).[UniProtKB/Swiss-Prot Function]