

## Product datasheet for MR222968L3

### Ctnnd2 (NM\_008729) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ctnnd2 (NM_008729) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Ctnnd2
Synonyms:	Catnd2; neurojugin; Nprap
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222968).
Restriction Sites:	SgfI-RsrII
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_008729
ORF Size:	3738 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_008729.2</a>
<b>RefSeq Size:</b>	5959 bp
<b>RefSeq ORF:</b>	3741 bp
<b>Locus ID:</b>	18163
<b>UniProt ID:</b>	<a href="#">O35927</a>
<b>Cytogenetics:</b>	15 B2
<b>Gene Summary:</b>	Has a critical role in neuronal development, particularly in the formation and/or maintenance of dendritic spines and synapses (PubMed:17993462) (PubMed:25807484). Involved in the regulation of canonical Wnt signaling (By similarity). It probably acts on beta-catenin turnover, facilitating beta-catenin interaction with GSK3B, phosphorylation, ubiquitination and degradation (PubMed:20623542). May be involved in neuronal cell adhesion and tissue morphogenesis and integrity by regulating adhesion molecules. Functions as a transcriptional activator when bound to ZBTB33 (PubMed:15282317).[UniProtKB/Swiss-Prot Function]

## Product images:



Circular map for MR222968L3