

Product datasheet for **MR222971L3V**

Fgf6 (NM_010204) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgf6 (NM_010204) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf6
Synonyms:	Fgf-6; HSTF-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010204
ORF Size:	624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR222971).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010204.1 , NP_034334.1
RefSeq Size:	1232 bp
RefSeq ORF:	627 bp
Locus ID:	14177
UniProt ID:	P21658
Cytogenetics:	6 61.92 cM
Gene Summary:	Plays an important role in the regulation of cell proliferation, cell differentiation, angiogenesis and myogenesis, and is required for normal muscle regeneration.[UniProtKB/Swiss-Prot Function]



[View online »](#)