

Product datasheet for **MR222975L4V**

II31 (NM_029594) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | II31 (NM_029594) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | II31 |
| Synonyms: | 1700013B14Rik |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_029594 |
| ORF Size: | 489 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR222975). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_029594.1 |
| RefSeq Size: | 810 bp |
| RefSeq ORF: | 492 bp |
| Locus ID: | 76399 |
| UniProt ID: | Q6EAL8 |
| Cytogenetics: | 5 F |



[View online »](#)

Gene Summary:

Activates STAT3 and possibly STAT1 and STAT5 through the IL31 heterodimeric receptor composed of IL31RA and OSMR (PubMed:15184896). May function in skin immunity (PubMed:15184896). Enhances myeloid progenitor cell survival in vitro (PubMed:17379091). Induces RETNLA and serum amyloid A protein expression in macrophages (PubMed:25847241).[UniProtKB/Swiss-Prot Function]