

Product datasheet for **MR223002L4V**

Atp4b (NM_009724) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp4b (NM_009724) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atp4b
Synonyms:	AV080843
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009724
ORF Size:	885 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223002).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009724.2 , NP_033854.1
RefSeq Size:	1327 bp
RefSeq ORF:	885 bp
Locus ID:	11945
UniProt ID:	P50992
Cytogenetics:	8 5.89 cM
Gene Summary:	Required for stabilization and maturation of the catalytic proton pump alpha subunit and may also involved in cell adhesion and establishing epithelial cell polarity.[UniProtKB/Swiss-Prot Function]



[View online »](#)