

Product datasheet for **MR223007L3V**

Cacnb1 (NM_001159319) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Cacnb1 (NM_001159319) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Cacnb1 |
| Synonyms: | CAB1; Cchb1; Cchl1b1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001159319 |
| ORF Size: | 1437 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR223007). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001159319.1 , NP_001152791.1 |
| RefSeq Size: | 1755 bp |
| RefSeq ORF: | 1440 bp |
| Locus ID: | 12295 |
| UniProt ID: | Q8R3Z5 |
| Cytogenetics: | 11 61.5 cM |



[View online »](#)

Gene Summary:

Regulatory subunit of L-type calcium channels. Regulates the activity of L-type calcium channels that contain CACNA1A as pore-forming subunit (By similarity). Regulates the activity of L-type calcium channels that contain CACNA1C as pore-forming subunit and increases the presence of the channel complex at the cell membrane. Required for functional expression L-type calcium channels that contain CACNA1D as pore-forming subunit. Regulates the activity of L-type calcium channels that contain CACNA1B as pore-forming subunit (By similarity). [UniProtKB/Swiss-Prot Function]