

Product datasheet for **MR223013L3V**

Sgcg (NM_011892) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sgcg (NM_011892) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sgcg
Synonyms:	35kDa; 5430420E18Rik; A1642964; gamma-SG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011892
ORF Size:	873 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223013).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011892.3 , NP_036022.1
RefSeq Size:	3515 bp
RefSeq ORF:	876 bp
Locus ID:	24053
UniProt ID:	P82348
Cytogenetics:	14 D1
Gene Summary:	Component of the sarcoglycan complex, a subcomplex of the dystrophin-glycoprotein complex which forms a link between the F-actin cytoskeleton and the extracellular matrix. [UniProtKB/Swiss-Prot Function]



[View online »](#)