

Product datasheet for **MR223020L4V**

C330027C09Rik (NM_172616) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C330027C09Rik (NM_172616) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	C330027C09Rik
Synonyms:	AA408511; AU018569; Cip2a; Kiaa1524
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172616
ORF Size:	2721 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223020).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172616.2 , NP_766204.2
RefSeq Size:	3974 bp
RefSeq ORF:	2724 bp
Locus ID:	224171
UniProt ID:	Q8BWY9
Cytogenetics:	16 B5
Gene Summary:	Promotes anchorage-independent cell growth and tumor formation.[UniProtKB/Swiss-Prot Function]



[View online »](#)