

## Product datasheet for **MR223021L4V**

### Ets1 (NM\_001038642) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ets1 (NM_001038642) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ets1
Synonyms:	AI196000; AI448617; D230050P06; Ets-1; Tpl1; vs
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001038642
ORF Size:	1059 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001038642.1</a>
RefSeq Size:	4799 bp
RefSeq ORF:	1062 bp
Locus ID:	23871
Cytogenetics:	9 17.97 cM



[View online »](#)

**Gene Summary:**

Transcription factor. Directly controls the expression of cytokine and chemokine genes in a wide variety of different cellular contexts. May control the differentiation, survival and proliferation of lymphoid cells. May also regulate angiogenesis through regulation of expression of genes controlling endothelial cell migration and invasion (By similarity). [UniProtKB/Swiss-Prot Function]