

Product datasheet for **MR223023L2V**

Braf (NM_139294) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Braf (NM_139294) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Braf
Synonyms:	9930012E13Rik; AA120551; AA387315; AA473386; B-raf; Braf-2; Braf2; C87398; C230098H17
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_139294
ORF Size:	2412 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223023).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139294.5 , NP_647455.3
RefSeq Size:	9728 bp
RefSeq ORF:	2415 bp
Locus ID:	109880
UniProt ID:	P28028
Cytogenetics:	6 18.43 cM
Gene Summary:	Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)