

## Product datasheet for **MR223029L3V**

### Crx (NM\_001113330) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Crx (NM_001113330) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Crx
Synonyms:	Crx1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113330
ORF Size:	969 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223029).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001113330.1</a> , <a href="#">NP_001106801.1</a>
RefSeq Size:	2903 bp
RefSeq ORF:	972 bp
Locus ID:	12951
UniProt ID:	<a href="#">O54751</a>
Cytogenetics:	7 8.6 cM



[View online »](#)

**Gene Summary:**

Transcription factor that binds and transactivates the sequence 5'-TAATC[CA]-3' which is found upstream of several photoreceptor-specific genes, including the opsin genes. Acts synergistically with other transcription factors, such as NRL, RORB and RAX, to regulate photoreceptor cell-specific gene transcription. Essential for the maintenance of mammalian photoreceptors.[UniProtKB/Swiss-Prot Function]