

## Product datasheet for **MR223054L3V**

### Ctbp1 (NM\_013502) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctbp1 (NM_013502) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctbp1
Synonyms:	BARS; CtBP3/BARS; D4S115h; D5H4S115; D5H4S115E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013502
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223054).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013502.3</a> , <a href="#">NP_038530.1</a>
RefSeq Size:	2305 bp
RefSeq ORF:	1326 bp
Locus ID:	13016
UniProt ID:	<a href="#">O88712</a>
Cytogenetics:	5 17.52 cM



[View online »](#)

**Gene Summary:**

Corepressor targeting diverse transcription regulators such as GLIS2 or BCL6. Has dehydrogenase activity. Involved in controlling the equilibrium between tubular and stacked structures in the Golgi complex. Functions in brown adipose tissue (BAT) differentiation. [UniProtKB/Swiss-Prot Function]