

Product datasheet for **MR223067L4V**

Ednrb (NM_007904) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ednrb (NM_007904) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ednrb
Synonyms:	ET-B; ET-BR; ETb; ETR-; ETR-b; Sox10; Sox10m1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007904
ORF Size:	1326 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223067).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007904.3 , NP_031930.1
RefSeq Size:	3989 bp
RefSeq ORF:	1329 bp
Locus ID:	13618
UniProt ID:	P48302
Cytogenetics:	14 53.05 cM



[View online »](#)

Gene Summary:

This gene encodes a member of the G-protein coupled receptor family. It encodes a receptor for endothelins, peptides that are involved in vasoconstriction. The encoded protein activates a phosphatidylinositol-calcium second messenger system and is required for the development of enteric neurons and melanocytes. Gene disruption causes pigmentation anomalies, deafness, and abnormal dilation of the colon due to defects of neural crest-derived cells. Mutations in this gene are found in the piebald mouse, and mouse models of Hirschsprung's disease and Waardenburg syndrome type 4. Renal collecting duct-specific gene deletion causes sodium retention and hypertension. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013]