

Product datasheet for **MR223078L4V**

Slc6a8 (NM_001142809) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc6a8 (NM_001142809) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc6a8
Synonyms:	AA589632; Creat; CRT; CRTR; CT1; CTR5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001142809
ORF Size:	1905 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223078).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001142809.1 , NP_001136281.1
RefSeq Size:	3977 bp
RefSeq ORF:	1908 bp
Locus ID:	102857
UniProt ID:	Q8VBW1
Cytogenetics:	X 37.38 cM
Gene Summary:	Required for the uptake of creatine. Plays an important role in supplying creatine to the brain via the blood-brain barrier.[UniProtKB/Swiss-Prot Function]



[View online »](#)