

## Product datasheet for MR223081L4

## Bcl2l11 (NM\_009754) Mouse Tagged Lenti ORF Clone

## **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** Bcl2l11 (NM 009754) Mouse Tagged Lenti ORF Clone

Tag: mGFP Symbol: Bcl2l11

Synonyms: 1500006F24Rik; bcl2-L-11; Bim; Bod

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

Sgfl-Mlul


E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(MR223081).

Sequence:

Restriction Sites: Cloning Scheme:





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_009754

ORF Size: 330 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 009754.3</u>, <u>NP 033884.1</u>

 RefSeq Size:
 4782 bp

 RefSeq ORF:
 333 bp

 Locus ID:
 12125


 UniProt ID:
 054918

 Cytogenetics:
 2 F1

**Gene Summary:** Induces apoptosis and anoikis. The isoforms vary in cytotoxicity with isoform BimS being the

most potent and isoform BimEL being the least potent.[UniProtKB/Swiss-Prot Function]

## **Product images:**



Circular map for MR223081L4