

## Product datasheet for **MR223094L3V**

### Clec9a (NM\_172732) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Clec9a (NM_172732) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Clec9a
Synonyms:	9830005G06Rik; DNGR-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172732
ORF Size:	714 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223094).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172732.3</a> , <a href="#">NP_766320.2</a>
RefSeq Size:	1156 bp
RefSeq ORF:	717 bp
Locus ID:	232414
UniProt ID:	<a href="#">Q8BRU4</a>
Cytogenetics:	6 F3



[View online »](#)

**Gene Summary:**

Functions as an endocytic receptor on a small subset of myeloid cells specialized for the uptake and processing of material from dead cells. Recognizes filamentous form of actin in association with particular actin-binding domains of cytoskeletal proteins, including spectrin, exposed when cell membranes are damaged, and mediate the cross-presentation of dead-cell associated antigens in a Syk-dependent manner.[UniProtKB/Swiss-Prot Function]