

Product datasheet for **MR223110**

Wdr70 (NM_001081402) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wdr70 (NM_001081402) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wdr70
Synonyms:	4833422F06Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR223110 representing NM_001081402
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGCACTCGGAACAGCAAGTGACAGGCGCAGACACGGCGGGACCCGACCCTCAGCTGGCGGTCA
CCATGGGCTTCACAGGATTCGGGAAGAAAGCTCGCACATTTGACTTGGAGGCCATGTTTCGAGCAAACCTG
AAGGACAGCGGTGAAAGAAGCCGAAAAACTGGAAGCAAGAGAAAAAGAGGAAGAAATGAACAGAGAA
AAAGAATTAAGAAAACAAATTGAAGATATGGAGCCAGCACCGTCAAGCTCGAGTGCAGCCAGAGAGCGCT
CACAGTCAATGCAGGGACACAAGCAGCAGTGACAGTGAGAGTGACGACAGTTCTGATTCATCTGACGA
TGAGCTAATCGGCCCTCCTTTACCCCTAAAATGGTAGGAGAATCCGTAACACTGTGGACGAAGGCACC
CTTGGCCCTTTACCTCCACCTCTGTGTGAAGAAGGAGAAGATGATGATGACGAGTTAGACGATGAAG
GAGAAGAGGATAACCCTGTCCAGAGGATTCCTGACTCGCATGAGATAAACTGAGGCATGGCACTAAAAC
AGTTTCTGCTCTGGGTCTGGACCCTCAGGTGCCCGTTTGGTCACTGGAGGATATGATTACGATGTTAAA
TTTTGGGATTTTCTGGAATGGATGCTTCTTCAAGGCATTTTCGATCCCTTCAACCTTGTGAATGCCATC
AGATCAAGTCGTTACAGTACAGTAACACAGGAGACATGATTCTTGTGTTTCTGGAAGCTCACAGGCCAA
GGTTATTGACAGAGATGGTTTTGAAGTCATGGAGTGTATAAAGGGAGACCAGTATATCGTAGACATGGCC
AACACCAAGGGCCACACAGCAATGCTTCATACTGGCTCATGGCATCCAAAAATAAAGGAGAATTTATGA
CTTGCTCAAATGATGCTACTGTGAGGCTCTGGGAAGTTGAAAAATCCAAAGAAGCAAAAAAGTGTTTTTAA
ACCTCGGACAATGCAGGGTAAAAAGTCAATCCACGACGTGCACGTACAGCAGAGACGGGAACCTGTG
GCCGCTGCCTGCCAGAACGGAAGCATTCCAGATCTGGGACCGGAACCTGACAGTTCACCCTAAATCCACT
ACAAGCAGGCTCACGACCCTGGCACAGACACGTCTTGTGTGGCTTTTTCCCTATGATGGTAATGCCTTGC
CTCTCGTGGAGGCGATGATACATTAATAATATGGGATGTCCGGCAATTAATAAGCCATTATTTTCAGCC
TCAGACCTTCTACCTTGTTCCTCAATGACCGACTGCTGCTTACGCTTACGACAAACTCATAGTTACCG
GTACCTCTGTCCAAAGGGGTGTGGCAGTGGCAAATCTGTCTTCTGGAGCGCAGGACTTTCCAAAGGGT
CTATGAAATCCACATCACAGATGCGAGTGTGGTCCGCTGCCTGTGGCATCCAAAGCTGAACCAGATCATG
GTTGGAATGGAATGGACTGGCAAAAGTGTACTACGACCCCAACAAGAGTCAGAGGGGAGCAAACTAT
GTGTGGTAAAACCCAGAGGAAGGCAAAACAAGCTGAGACACTAACCCAGGACTACATCATACCCCTCA
TGCCTTGCCCATGTTCCGTGAGCCCGCCAGCGAAGTACGAGGAAACAGCTGGAGAAGGACAGACTGGAC
CCCCTGAAGTCCACAAACCAGAGCCTCCTGTAGCAGGCCAGGTCGTGGTGGCCGAGTTGGAACCCACG
GGGGCACCTGTCTCTACATTGTGAAGAACATCGCTTTGGACAAGACCAGCAGACGAATCCCGGGGA
AGCCATTTGCGCCATGCCAAAGCAGCAGAAGATAACCCATATTGGGTTTCGCCAGCATATTCAAAAC
CAGCAAAAACCATGTTTGCCCAAGTGGAGTCTGATGATGAAGAGTCAAAGAATGAACCAGAATGGAAAA
AACGTAATAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223110 representing NM_001081402
Red=Cloning site Green=Tags(s)


```
MEHSGTSEVTGADTAGPDPQLAVTMGFTGFGKARTFDLEAMFEQTRRTAVERSRLTLEAREKEEEMNRE
KELRKQIEDMEPAPSSSSAARERSQSSCRDTSSSDSEDDSSDSDDELIGPPLPPKVMGESVTTVDEGT
LGPLPPPLCEEGEDDDDELDEGEEDNPVQRIPDSHEITLRHGKTKVSALGLDPSGARLVTGGYDYDVK
FWDFAGMDASFKAFRSLQPCECHQIKSLQYSNTGDMILVVSQSSQAKVIDRDGFVMECIKGDQYIVDMA
NTKGHTAMLHTGSWHPKIKGEFMTCSNDATVRLWEVENPKKQKSVFKPRMQGKVIPTTCTYSRDGNLV
AAACQNGSIQIWRNLTVHPKFHYKQAHDPGTDTSVAFSYDGNVLASRGDDTLKLWDVRFNKPLFSA
SDLPTLFPMTDCCFSPDDKLIIVTGTSVQRGCGSKLVFLERRTFQRVYIEIHTDASVVRCLWHPKLNQIM
VGTGNGLAKVYYPNKSQRGAKLCVVKTRKAKQAETLTQDYIITPHALPMFREPRQRSTRKQLEKDRLD
PLKSHKPEPPVAGPGRGVRVGTGGTSSYIVKNIALDKTDDSNPREAILRHAKAAEDNPYWSPAYSKT
QPKTMFAQVESDDEESKNEPEWKKRKI
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9045_f11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001081402

ORF Size: 1971 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001081402.1](#), [NP_001074871.1](#)

RefSeq Size: 2151 bp

RefSeq ORF: 1974 bp


Locus ID: 545085

UniProt ID: [Q3TWF6](#)

Cytogenetics: 15 A1

MW: 73.5 kDa

Product images:



Circular map for MR223110