

Product datasheet for **MR223112L4V**

Cdh7 (NM_172853) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdh7 (NM_172853) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdh7
Synonyms:	9330156F07Rik; CDH7L; CDH7L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172853
ORF Size:	2355 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223112).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172853.2 , NP_766441.1
RefSeq Size:	3353 bp
RefSeq ORF:	2358 bp
Locus ID:	241201
UniProt ID:	Q8BM92
Cytogenetics:	1 E2.1



[View online »](#)

Gene Summary:

This gene encodes a member of the cadherin family of calcium-dependent glycoproteins that mediate cell adhesion and regulate many morphogenetic events during development. The encoded preproprotein is further processed to generate a mature protein. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene is located adjacent to a related cadherin gene on chromosome 1. [provided by RefSeq, Oct 2015]