

## Product datasheet for **MR223148L3V**

### **Gpm6b (NM\_001177959) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Gpm6b (NM_001177959) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gpm6b
Synonyms:	AI593561; Gpm6; M6B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177959
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223148).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001177959.1</a> , <a href="#">NP_001171430.1</a>
RefSeq Size:	4828 bp
RefSeq ORF:	798 bp
Locus ID:	14758
UniProt ID:	<a href="#">P35803</a>
Cytogenetics:	X F5



[View online »](#)

**Gene Summary:**

May be involved in neural development. Involved in regulation of osteoblast function and bone formation. Involved in matrix vesicle release by osteoblasts; this function seems to involve maintenance of the actin cytoskeleton. May be involved in cellular trafficking of SERT and thereby in regulation of serotonin uptake.[UniProtKB/Swiss-Prot Function]