

## Product datasheet for **MR223154L3V**

### Tnnt2 (NM\_001130176) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnnt2 (NM_001130176) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnnt2
Synonyms:	cTnT; Tnt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001130176
ORF Size:	903 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223154).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001130176.1</a> , <a href="#">NP_001123648.1</a>
RefSeq Size:	1150 bp
RefSeq ORF:	906 bp
Locus ID:	21956
UniProt ID:	<a href="#">P50752</a>
Cytogenetics:	1 59.32 cM
Gene Summary:	Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity. [UniProtKB/Swiss-Prot Function]



[View online »](#)