

## Product datasheet for **MR223171L4V**

### PPP4R3A (NM\_001160214) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppp4r3a (NM_001160214) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ppp4r3a
Synonyms:	1110034C04Rik; BC064465; mKIAA2010; Smek1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001160214
ORF Size:	2499 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223171).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001160214.1</a> , <a href="#">NP_001153686.1</a>
RefSeq Size:	4594 bp
RefSeq ORF:	2502 bp
Locus ID:	68734
UniProt ID:	<a href="#">Q6P2K6</a>
Cytogenetics:	12 E


[View online »](#)

**Gene Summary:**

Regulatory subunit of serine/threonine-protein phosphatase 4. May regulate the activity of PPP4C at centrosomal microtubule organizing centers. The PPP4C-PPP4R2-PPP4R3A PP4 complex specifically dephosphorylates H2AFX phosphorylated on 'Ser-140' (gamma-H2AFX) generated during DNA replication and required for DNA DSB repair (By similarity).  
[UniProtKB/Swiss-Prot Function]