

Product datasheet for MR223174

Nfasc (NM_001160318) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfasc (NM_001160318) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfasc
Synonyms:	AA387016; D430023G06Rik; mKIAA0756; NF
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR223174 representing NM_001160318 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCAGGCAGCAGGCGCCACCCTGGGTCCACATAGCCCTCATCTCTTCTCCTCAGCCTCGGAGGGG
CCATCGAGATCCGATGGACCTGACCCAACCCCAACTATCACCAAGCAGTCCGTGAAGGACCACATCGT
GGACCCTCGAGATAACATCCTGATTGAATGTGAAGCTAAAGGCAACCCCGCCCCAGTTTTCTACTGGACT
CGCAACAGCAGATTCTTCAACATTGCCAAGGACCCACGGGTGTCCATGAGGAGGAGATCTGGACCTTGG
TGATCGACTTCCGAGTGGTGGGCGGCCTGAGGAATACGAAGGGGAGTACCAGTGCTTGGCCGGAACAA
ATTTGGCACTGCACCTTAGCAACCGCATCCGCTGCAGGTGCCAAATCTCCCTGTGGCCCAAGGAAAAAC
CTAGACCCCGTCGTGGTTCAAGAGGGTGGCCCTTGACACTACAGTGCAACCCCCACCCGGCCTCCCGT
CCCCCGTCATCTTCTGGATGAGCAGTCCATGGAGCCATCACCAGGACAAGCGTGTCTCCAGGGTCA
CAACGGGGACCTGTACTTCTCAACGTCATGCTGCAGGACATGCAGACCGACTACAGTGCACCGCGCGC
TTTCACTTCACCCACACCATTTCAGCAGAAGAACCCTTCACCTCAAGTCTCACCACCCGAGGAGTTG
CAGAAAGAAGCCAGCTTCATGTATCCCCAGGGCACATCGAGCAGTCAAGTGGTCTCCGTGGCATGGA
CCTGTCTGGAATGCATTGCCTCTGGCGTCCCAACACCAGACATTGCATGGTACAAGAAAGGTGGGGAC
CTCCCATCTAACCAAGGCCAAGTTCGAGAAGTTAATAAGGCTCTGCGCATCACCATGTCTCTGAAGAGG
ACTCTGGGAGTATTTCTGCCTGGCCTCCAACAAGATGGGCAGCATCCGGCACACGATCTCGGTGAGAGT
AAAGGCTGCTCCATACTGGCTGGATGAGCCCAAGAACCCTGATCCTGGCTCCTGGGGAAGATGGGAGGCTG
GTATGCCGAGCCAATGGGAACCCGAAGCCGACCGTGCAGTGGTGAATGGAGAGCCTTTACAATCGG
CACCACCAATCCAACCGTGAGGTAGGTGAGACACTATCATCTCCGGGATACTCAGATCAGCAGCAG
GGCAGTGTACCAATGTAATACATCCAATGAACATGGCTACCTGCTGGCCAATGCCTTCGTCAGCGTGT
GATGTACCCCTCGGATGCTGTCTGCCCAACCAGCTCATCAGGGTATCCTTTATAACCGGACACGGC
TGGACTGTCGGTCTTTGGGTCTCCCATCCAACACTCCGATGGTTAAGAATGGCAAGGAAGCAACCT
GGATGGCGGTAACCTACCAGTCTACGAAAACGGCAGTCTAGAAATCAAGATGATTCGCAAGAGGACCAA
GGCATCTACACCTGTGTGGCCACCAACATCCTGGGCAAAGCCGAAATCAAGTCCGCTGGAGGTCAAAG



[View online »](#)

ACCCACCCAGGATCTACAGGATGCCCGAGGACCAGGTGGCCAAGAGGGGCACCACGGTGCAGCTGGAGTG
CCGCGTGAACATGACCCTCCTTGAAGCTCACAGTCTCCTGGCTGAAGGACGATGAGCCACTCTACATT
GGAAACAGGATGAAGAAGGAAGATGACTCCCTGACGATCTTCGGAGTGGCAGAGCGGGACCAGGGCAGTT
ACACGTGTATGGCCAGCACCCGAGCTGGACCAGGACCTGGCAAAGGCTACCTCACTGTTCTAGGGCGACC
AGACCCGACCCAGGACCTGGAGCTCACTGACCTGGCTGAGAGGAGTGTGAGGCTGACCTGGATCCAGGG
GATGACAACAACAGCCCTATCACAGACTACGTCGTTCAGTTTGAAGAGGACCAGTTCCAACCAGGGGTGT
GGCATGACCACCTCAGGTTCCAGGCAGCGTCAACTCAGCCGCTCCATCTGTCCCATATGTCAACTA
CCAATTCAGAGTCATCGCTGTCAACGAGGTTGGGAGCAGCCACCCAGCCTTCCATCCGAGCGGTACCCGA
ACCAGTGGGGCACCCCTGAATCTAATCCAGTGTGAAGGGCGAAGGGACAAGAAAGAACAATATGG
AGATCACGTGGACGCCTATGAATGCTACCTCTGCCTTTGGCCCCAACCTACGCTACATTGTCAAGTGGCG
ACGGAGAGAGACCCGAGAGACTTGAACAATGTACAGTGTGGGCTCTCGCTACGTGGTGGGCAGACG
CCTGTCTACGTTCCCTATGAGATCCGAGTCCAGGCTGAAAATGACTTTGGGAAAGCCCCGAGCCTGACA
CCATCATTGGGTACTCCGGAGAAGATTATCCAGGGCTGCGCCCACTGAAGTTAAATCCGAGTCTGAA
CAGCACAGCCATCAGCCTTCAGTGAACCGAGTCTACTCTGACACGGTCCAGGGCCAGCTCAGAGAGTAT
CGAGCTTACTACTGGAGGAAAGCAGTTTGTGAAGAACCTGTGGGTGTCTCAGAAGAGACAGCAGGCCA
GCTTCCCTGGTGACCGTCCCCGGGGCGTGGTGGCCCGCTGTTCCCTACAGTAACTACAAGCTGGAGAT
GGTGGTGGTCAATGGGAGAGGTGACGGGCTCGAAGTGAACCAAGGAATTCACCACCCAGAAAGGAGTA
CCCAGTGGCCCCAGGCGTTTCAAGTCCGACAGCCCACTGGAGACCATCAACCTGGAGTGGGACCACC
CAGAGCACCCCAACGGAATCCTGATTGGATACATCCTCAGATACGTGCCCTTTAATGGAACCAACTGGG
AAAGCAGATGGTGGAAAACCTTCTCTCCAATCAGACCAAGTTCTCTGTGCAGAGAGCAGACCCAGTGTG
CGTTACCGCTTCTCCCTCAGTCCAGGACACAGGTGGGCTCTGGAGAAGCAGCCACAGAGGAGTCCCCAG
CACCTCCAAATGAAGCTACTCCAACGACGTTACACCAATAACCAGGCAGACATCGCCACCCAGGGCTG
GTTTCATCGGGCTCATGTGTGCCATTGCCCTTCTGGTGTGATCCTTCTCATCGTCTGCTTTCATCAAGAGG
AGTCGAGGTGGCAAGTACCCAGTCCGAGGAAAAGAAGGATGTCCCTTGGTCTGAAGACCCAAAGAAG
AAGATGGCTCATTGACTACAGTGTGAGGACAACAAGCCCTGCAGGGCAGCCAGACATCTCTGGATGG
CACCATCAAGCAGCAGGAGAGCGATGACAGCCTGGTGGACTATGGCGAAGGCGGGCAGGGCCAGTTCAAT
GAAGATGGCTCCTTTATTGGCCAGTACTGTCAAAAAGGACAAGGAGGAAACGGAGGGCAATGAGAGCT
CAGAGGCCACATCACCAGTCAATGCCATCTATTCCCTTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223174 representing NM_001160318
 Red=Cloning site Green=Tags(s)

MARQQAPPWVHIALILFLLSLGGAIEIPMDLTQPPTITKQSVKDHI V D P R D N I L I E C E A K G N P A P S F H W T
 R N S R F F N I A K D P R V S M R R R S G T L V I D F R S G G R P E E Y E G E Y Q C F A R N K F G T A L S N R I R L Q V S K S P L W P K E N
 L D P V V V Q E G A P L T L Q C N P P P G L P S P V I F W M S S S M E P I T Q D K R V S Q G H N G D L Y F S N V M L Q D M Q T D Y S C N A R
 F H F T H T I Q Q K N P F T L K V L T T R G V A E R T P S F M Y P Q G T S S S Q M V L R G M D L L L E C I A S G V P T P D I A W Y K K G G D
 L P S N K A K F E N F N K A L R I T N V S E E D S G E Y F C L A S N K M G S I R H T I S V R V K A A P Y W L D E P K N L I L A P G E D G R L
 V C R A N G N P K P T V Q W M V N G E P L Q S A P P N P N R E V A G D T I I F R D T Q I S S R A V Y Q C N T S N E H G Y L L A N A F V S V L
 D V P P R M L S A R N Q L I R V I L Y N R T R L D C P F F G S P I P T L R W F K N G Q G S N L D G G N Y H V Y E N G S L E I K M I R K E D Q
 G I Y T C V A T N I L G K A E N Q V R L E V K D P T R I Y R M P E D Q V A K R G T T V Q L E C R V K H D P S L K L T V S W L K D D E P L Y I
 G N R M K K E D D S L T I F G V A E R D Q G S Y T C M A S T E L D Q D L A K A Y L T V L G R P D R P R D L E L T D L A E R S V R L T W I P G
 D D N N S P I T D Y V V Q F E E D Q F Q P G V W H D S R F P G S V N S A V L H L S P Y V N Y Q F R V I A V N E V G S S H P S L P S E R Y R
 T S G A P P E S N P S D V K G E G T R K N N M E I T W T P M N A T S A F G P N L R Y I V K W R R R E T R E T W N N V T V W G S R Y V V G Q T
 P V Y V P Y E I R V Q A E N D F G K G P E P D T I I G Y S G E D Y P R A A P T E V K I R V L N S T A I S L Q W N R V Y S D T V Q G Q L R E Y
 R A Y Y W R E S S L L K N L W V S Q K R Q Q A S F P G D R P R G V A R L F P Y S N Y K L E M V V N G R G D G P R S E T K E F T T P E G V
 P S A P R R F R V R Q P N L E T I N L E W D H P E H P N G I L I G Y I L R Y V P F N G T K L G K Q M V E N F S P N Q T K F S V Q R A D P V S
 R Y R F S L S A R T Q V G S G E A A T E E S P A P P N E A T P T A A Y T N N Q A D I A T Q G W F I G L M C A I A L L V L I L L I V C F I K R
 S R G G K Y P V R E K K D V P L G P E D P K E E D G S F D Y S D E D N K P L Q G S Q T S L D G T I K Q Q E S D D S L V D Y G E G G E G Q F N
 E D G S F I G Q Y T V K K D K E E T E G N E S S E A T S P V N A I Y S L A


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001160318

ORF Size: 3471 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001160318.1](#), [NP_001153790.1](#)

RefSeq Size: 9647 bp

RefSeq ORF: 3474 bp

Locus ID: 269116

Cytogenetics: 1 57.42 cM

MW: 130.7 kDa

Gene Summary:

This gene encodes an L1 family immunoglobulin cell adhesion molecule with multiple IGcam and fibronectin domains. The protein functions in neurite outgrowth, neurite fasciculation, and organization of the axon initial segment (AIS) and nodes of Ranvier on axons during early development. Both the AIS and nodes of Ranvier contain high densities of voltage-gated Na⁺ (Nav) channels which are clustered by interactions with cytoskeletal and scaffolding proteins including this protein, gliomedin, ankyrin 3 (ankyrin-G), and betaIV spectrin. This protein links the AIS extracellular matrix to the intracellular cytoskeleton. This gene undergoes extensive alternative splicing, and the full-length nature of some variants has not been determined. [provided by RefSeq, May 2009]