

Product datasheet for **MR223210L3V**

Tomm40 (NM_001109748) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tomm40 (NM_001109748) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tomm40
Synonyms:	AW539759; Mom35; Tom40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001109748
ORF Size:	1086 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223210).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001109748.1 , NP_001103218.1
RefSeq Size:	1669 bp
RefSeq ORF:	1086 bp
Locus ID:	53333
UniProt ID:	Q9QYA2
Cytogenetics:	7 A3
Gene Summary:	Channel-forming protein essential for import of protein precursors into mitochondria. [UniProtKB/Swiss-Prot Function]



[View online »](#)