

Product datasheet for **MR223213A1V**

Mouse Bmp3 (NM_173404) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Bmp3 (NM_173404) AAV Particle
Tag:	Myc-DDK
Symbol:	Bmp3
Synonyms:	9530029I04Rik; D630004D15R
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_173404.3, NP_775580.1</u>
RefSeq Size:	3112 bp
RefSeq ORF:	1407 bp
Locus ID:	110075
UniProt ID:	<u>Q8BHE5</u>
Cytogenetics:	5 48.24 cM
MW:	52.8 kDa