

Product datasheet for **MR223266**

Trim25 (NM_009546) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trim25 (NM_009546) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Trim25
Synonyms:	AA960166; AL022677; EFP; Zfp147
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR223266 representing NM_009546
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGGAGCTGAATCCTCTGGCCGAGGAGCTGTCTGCTCCGTCTGCCTGGAGCTCTTCAAGAACCTG
 TCACCACTCCGTGCGGCCACAATTTCTGCATGTATGCCTGGATGAGACGTGGGTCTCCAGGGCCCGCC
 GTACCGATGCCCGCAATGTCGCAAAGTGTACCAGGTGCGCCCGCAGCTGCAAAAGAACACGGTAATGTGC
 GCGGTGGTGGAGCAGTTCTGCAAGGCGAGCAGGCGCGGACCCCGGTGGACGACTGGACGCCTCTGCC
 GCTTCTCCGCCTCCAGCGCAGCCACCCAGGTGGCCTGCGACCACTGCCTGACGGAGATCGCAGTGAAGAC
 GTGCCTCGTGTGCATGGCCTCTTTCTGCCAGGAGCACTTGAGGCCTCACTTCGACAGTCCGGCCTCCAG
 GATCACCCCTACAGTCGCCATTCGTGACCTTCTGCGCCGCAATGTACCCAGCACAACCGCCTTCGGG
 AATTGTTTTGTCGGAGCATGGCGAGTGCATCTGCCATATCTGCTTGGTAGAGCACAAGACCTGCTCCCC
 TACGACCTAAGTCAAGCCAGCGCCGACCTGGAGTACAACTGAGGAATAAACTCACTATCATGCACAGC
 CACATCAATGGGGCAACCAAGCATTGGAGGATGTGAGATCCAAGCAACAGTGTGTGCAGGATTCTATGA
 AGAGGAAGATGGAACAGCTGAGACAGGAGTATATGGAAATGAAGGCTGTCATCGACGCAGCAGAGACCAG
 CTCCTTGCGAAGTTAAAGGAGGAGGAGAAGAGAGTCTACGGAAAATTCGACACCATCTACCAGGTTCTC
 GTCAAGAAGAAGAGTGAAGTGCAGAACTGAAGGCAGAGGTTGAGCTCATTATGGACAAAGGGGATGAGT
 TTGAGTTTCTGGAGAAAGCTGCAAAATGCAAGGAGAGTCAACAAAACCAGTCTACATCCCTAAGATAGA
 TCTGGATCATGATTTGATAATGGGAATCTACCAAGGGCGCGGACCTCAAGAGTGAAGTGAAGCACTCC
 ATTAAGAAGCTACAGAAGAAGTCTGAGGAGCACAATGGCTCAGGTAACAAGGGAGACCAGACACAGTCCA
 CTTTTAAACCCGTACAACCTTCGAAGAAGACCATTCAAGAAAAGAAGACCAAGAAAACCCCTGTTGCTCC
 CGGCCCTCCCTCCCACTTCTCACCTAACAGTTGCCACCTTTGGAGCTCCTGGCCAATCACTGGATTCA
 AAAGCAACTTCCCTGATGCCGCACCAAGCCAGTGCTGCACAGCCAGACTCTGTAGGAGTCAAGGCTA
 AGGTGCTAGAAAATCTTAACCAAGTCCAGAACGGAGCTTCTGGAGTATTTCTGAAAAGTCACTCTCGA
 CTACAATACCGCCACAACAAAGTGTCTGTCAAACAAGTACACCACCGCCTCTGTGTCCGATGGGCTC
 CAGCACTACCGGTCCCATCCCCAGAGGTTACCTACTGCTCTCAGGTCTGGGGCTGCACTGCTACAAGA
 ACGGCATTCACTACTGGGAGGTGGAACACAGAAGAACAATCTGCGGCGTAGGCATCTGCTACGGCAG
 CATGGAACGGCAGGGCCCGAGAGCCGGCTGGGCCGAACCCCAACTCATGGTGTGTGGAGTGGTTCAAT
 AACAAGATCTCTGCCTGGCACAACAACGTGGAGAAGACTCTGCCCTCCACCAAGGCCACTCGAGTCGGCG
 TGCTGCTCAACTGTGACCACGGCTTCGTCAATTTCTTTGCGGTACCCGAAAAGGTGCACCTGATGTACAA
 GTTCAAGGTGGACTTCACCGAAGCTCTGTACCCAGCCTTCTGGGTGTTCTCGGCCGTACCACCTCTCT
 ATCTGCTCCAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223266 representing NM_009546
 Red=Cloning site Green=Tags(s)

```


MAELNPLAEELSCSVCLELFKEPVTTPCGHNFMSCLDETWVVQGGPPYRCPQCRKVVYQVRPQLQKNTVMC
AVVEQFLQAEQARTPVDDWTPPARFSASSAATQVACDHCLTEIAVKTCLVCMASFQEHLPFHDFSPAFQ
DHPLQSPIRDLLRRKCTQHNRLRELFCEPHGECICHICLVHKTCSPPTLSQASADLEYKLRNKLTIMHS
HINGATKALEDVRSKQQCVQDSMKRKMQLRQEYMEMKAVIDAAETSSLRKLKEEEKRVYGFDTIYQVL
YKKKSEMQLKAEVELIMDKGDEFEFLEKAAKLQGESTKPVYIPKIDLDHDLIMGIYQGAADLKSELKHS
IKKLQKKSEEHNHSGNKGQDTQSTFKPVQPSKKTIQEKTKKTPVAGPPSPHFSPNKLPTFGAPGQSLDS
KATSPDAAPKASAAQPDVSGVKAKVLENFLTKSRTELLEYFVKVIFDYNTAHNKVSLSNKYTTASVSDGL
QHYRSHPRQRTYCSQVLGLHCYKNGIHYWEVELQKNNFCGVGICYGSMERQGPESRLGRNPNSWCVEWFN
NKISAWHNNVEKTL PSTKATRVGVLLNCDHGFVIFFAVTEKVHLMYKFKVDFTEALYPAFWVF SAGTTLS
ICSK
  
```


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_009546

ORF Size: 1902 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009546.2](#), [NP_033572.2](#)

RefSeq Size: 5596 bp

RefSeq ORF: 1905 bp

Locus ID: 217069

UniProt ID: [Q61510](#)

Cytogenetics: 11 C

MW: 72.2 kDa

Gene Summary: Functions as a ubiquitin E3 ligase and as an ISG15 E3 ligase. Involved in innate immune defense against viruses by mediating ubiquitination of DDX58 and IFIH1. Mediates 'Lys-63'-linked polyubiquitination of the DDX58 N-terminal CARD-like region which is crucial for triggering the cytosolic signal transduction that leads to the production of interferons in response to viral infection. Mediates 'Lys-63'-linked polyubiquitination of IFIH1. Promotes ISGylation of 14-3-3 sigma (SFN), an adapter protein implicated in the regulation of a large spectrum signaling pathway. Mediates estrogen action in various target organs. Mediates the ubiquitination and subsequent proteasomal degradation of ZFH3.[UniProtKB/Swiss-Prot Function]