

Product datasheet for MR223277L4

Rassf8 (NM_027760) Mouse Tagged Lenti ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Rassf8 (NM_027760) Mouse Tagged Lenti ORF Clone |
| Tag: | mGFP |
| Symbol: | Rassf8 |
| Synonyms: | 5133400D11Rik; AA815821; AW123240; Hoj1; mHoj-1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| E. coli Selection: | Chloramphenicol (34 ug/mL) |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR223277). |
| Restriction Sites: | SgfI-MluI |
| Cloning Scheme: | |

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

| | |
|-----------|-----------|
| ACCN: | NM_027760 |
| ORF Size: | 1260 bp |



[View online »](#)

