

Product datasheet for **MR223327**

Vmn1r177 (NM_206872) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Vmn1r177
Synonyms:	Vlrd12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR223327 representing NM_206872
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGTCTTCTCAGAATGGAGCCATGAAAACCACTGAGGAAGTGGCACTTCAGATTGTGTTGCTTTGTCAGT
TTGGGATAGGGACTATGGCTAACATTCTTCTGTTGTCCATAATTTCTCTCCACTTTTGACTGGATGTCA
GCTGAGACCCAAAGAAGTGATTTGCTAACTGCTGTGTCAGCATCTTAATTCTATTCATGACTGTG
TTTCGAAACAATGTGTTGTTGCTCAAAGGAATCCACTGAATGACCTTCAATGTAACTTGAGTACTTTG
TATACATGGTAGCAAGAAGCACAACTTGTGTTCCACCTGTGTCTTGAGCACCTATCAGTTTGTCACTCT
TGTTCTCGTAATTTCAATGAGCTGATGGTCAGAGGAAGAAGCCCCAAGGTCTAAGCTATTCCTGTTGT
AGTTGTTGGTTGTTGGTGTCTTATATAATGTCTACATTCCAATGAAAGTCAGTGGTCCACAGAACAGAC
ATAATAAAACCACTGCTAAAGTAAGTGGGCTGTTTTACATCTGGCTTCAGTGTGGCATTGTGTTCTT
GTCATTTGCCCATGATGCCATATTTATCAGCATCATGGTCTGGGCCAGTGTCTCCATGGTACTTCTCCTG
TATAGACACTACCAGAGACTGCAGTATATTTCCCTCCCCATCAAGGCAACAGAGGCTATGCTGAGATCA
GAGCAGCCCATAGCATTGTGATGCTTGTGGTCACATTTGTAAGCTTTTATCTTCTACACTGTATTTCTT
TATCTTACATTGTTTTGTAGACTCCCATCTGTGGTTGAGGCATATAGGTGAAGTTTAACTTCAAGC
TTCCCCACTATTTCTCCCTTGCTGTTGATCTTTAGGGATCCTAGTTATCATTGTTCTCTGATCTTCAACT
ATAAAAGCCAGTCATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR223327 representing NM_206872

Red=Cloning site Green=Tags(s)

MSSQNGAMKTTEEVALQIVLLCQFGIGTMANILLFVHNFSPLLTCQLRPKEVICTNMAVSSILILFMTV
FRNNVFVAQRNPLNDLQCKLEYFVYMVARSTNLCSTCVLSTYQFVTLVPGNFNELMVRGRSPKVISYSCC
SCWLFGLVYNVYIPMKVSGPQNRHNKTSTKSKWVCFTSGFSVGIVFLSFAHDAIFISIMVWASVSMVLLL
YRHYQRLQYIFPPHQGNRGYAEIRAAHSIVMLVVTFSFYLLHCISFILHIVFVDSHLWLRHIGEVLTSS
FPTISPLLLIFRDPSYHCSLIFNYKKPVI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mm9086_d07.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_206872

ORF Size:

927 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_206872.2, NP_996755.2

RefSeq Size: 930 bp


RefSeq ORF: 930 bp

Locus ID: 384572

Cytogenetics: 7 A3

MW: 35.2 kDa

Product images:



Circular map for MR223327