

Product datasheet for **MR22330L4V**

Gins1 (NM_027014) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gins1 (NM_027014) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gins1
Synonyms:	2810418N01Rik; Gins4; mKIAA0186; PSF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027014
ORF Size:	588 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223330).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027014.1 , NP_081290.1
RefSeq Size:	1100 bp
RefSeq ORF:	591 bp
Locus ID:	69270
UniProt ID:	Q9CZ15
Cytogenetics:	2 G3
Gene Summary:	Required for correct functioning of the GINS complex, a complex that plays an essential role in the initiation of DNA replication, and progression of DNA replication forks. GINS complex seems to bind preferentially to single-stranded DNA.[UniProtKB/Swiss-Prot Function]


[View online »](#)