

Product datasheet for **MR223332L3V**

Lrrc8b (NM_001033550) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Lrrc8b
Synonyms:	2210408K08Rik; mKIAA0231; R75581; Ta-Irrp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033550
ORF Size:	2409 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR223332).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001033550.2
RefSeq Size:	2841 bp
RefSeq ORF:	2412 bp
Locus ID:	433926
UniProt ID:	Q5DU41
Cytogenetics:	5 E5



Gene Summary:

Non-essential component of the volume-regulated anion channel (VRAC, also named VSOAC channel), an anion channel required to maintain a constant cell volume in response to extracellular or intracellular osmotic changes. The VRAC channel conducts iodide better than chloride and can also conduct organic osmolytes like taurine. Channel activity requires LRRC8A plus at least one other family member (LRRC8B, LRRC8C, LRRC8D or LRRC8E); channel characteristics depend on the precise subunit composition.[UniProtKB/Swiss-Prot Function]