

Product datasheet for **MR223336**

Ovch2 (NM_172908) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Ovch2
Synonyms:	9230106D23Rik; Ovtn
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR223336 representing NM_172908
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCCTATCAGCAAGGACAAGCTGATTTTAAATATTAGGAATGGTCTGCTTGAACAAGGTCATTCTGAAA
 CCCTTTCTCCATAAGAAATCCTGACTGTGGGCAGAGTTTGGTTAAGCCCCAACCTCAGAATTACTTTAG
 CCTTTTCAGTCGCATTGTTGGAGGAAGCCAGGTGGAAGGGTTCTATCCTTGGCAAGTGTCTCTGAAA
 CAAAAGCAGAAGCACATCTGTGGAGGGACCATCTCTTCACAGTGGGTATCACAGCCGCTCATTGCA
 TGGCTAACAGAAACATTGCATTAACCTTTGAATGTAAGTCTGGAGAATGACTTGAGCCAAGCCGAGCC
 AGGAGAGCAAAACATTGCCATTGAAACCATCATCATACCCGCGAGTTCTCCACCAGAAAACTATGATC
 TATGATATAGCCCTTCTGAAGATGGCTGGAACCTTCCAATTTGGACAGTTTGTGAGGCCTGTGTCTTC
 CAGAGCCAGGAGAACCTTTAATGCTGGATTTATTTGTACAACTGCAGGCTGGGTCGCTGTCCGAAGG
 TGGCAGGCTCCACAAGTCTTGACGAGGTAAACCTGCCAATTTTGAAGTCTGAGGAGGAATGTGAGGCCGT
 CTGTTGACCTTAAAGAATCCCATCTGGAAGACCTTCTCTGCACAGGATCCCTGATGGAGGGAGAG
 ATGCTTGTACAGGAGATTCAGGAGGGTCACTTATGTGCCAGAATAGGAAAGGGGCTGGACTCTGGCTGG
 AGTGACTTCTTGGGTTTGGGCTGTGGTCGAAGCTGGAGAAACAATGCAAGAAAGAAAGAGCAAGGATCT
 CCTGGAATCTTCACAGATCTTAGGAGAGTGCTTCCCTGGATCCTCAAACATATCCAGACTGGCCATCGAA
 GGAAGAGCACAAAGCTTTATGCAGTGAGCCTGATGGTCTCATCAGTGGGTCTGAGGGAGAGCTGCACTT
 CCCAGAAAGCCTCCACCTGTACTATGAGAGCAAGCAACTGTGTCTGGACCTTCTGGTTCAGAGGAC
 ATGCATATGTTGCTCAATTTGTACACTTGGATATAGAGTCTTGTACCCACAATTACCTGGCAATGTACT
 CTTTAGAAGACAGACTCGTTGGGAAATCTGTGGAGAAATCTCCCTTCATCCATTCTATTGGCTCCAG
 TTCCATCCGGCTGAGGTTTCATCTCTGATGCCACAGATTATGCCACTGGGTTAATCTCACCTATAAGCA
 CTTAAGCCGAGATACCATCCTGATTCGGCTGCAGGTCTCTAACTATCCTCTTTGAAGAAGGCACTATAC
 AGAGTCTTCACTATCCTGAAGACTACAGTAACATGGCTAGCTGTACCTGGATTTTCAAGCCCCAACTA
 TTGCCTAATTAAGCTGTATTTAGAGCCTAGAAATAGAAGAAATGGAGACTGCTCCTCTGACTATGTG
 ACTGTACACAGTGATGTCGAGAAGGAGAAGGAAATAGCTCGGTTCTGCGACTATGTGATACCCAGCCCTG
 TGCTGAGTTCCTCCAGTGATGCTCATCAGCTTCCAATCAGATGAAATGGCACTGCCAGGGGCTTCCA
 GGCTGACATTTCTTCTATTCTAGAGCAGATTTAAATATCTCCATATCAGAGGATGAATCAGTGCCTCTG
 GAGACATGGGACCTACCACCTGGAGCAATGAAATTTTCGATGCAGAGAGAGACACATACAAAACAC
 CATACGAAGAAGACATCGGGGAGATGCCAGCTATTGACAGTGGTCTTCTGAAGCAAGGGGAGAGGCGTGG
 AAAGCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223336 representing NM_172908
 Red=Cloning site Green=Tags(s)


MPISKDKLILILGMVLEQGHSETLSSIRNPDCGQSLVKPQPQNYFSLFSRIVGGSQVEKGSYPWQVSLK
 QKQKHICGGTIISQWVITAACHMANRNIALTLNVTAGEHDLQAEPGEQTLAIETIIHPQFSTRKPMI
 YDIALLKMAGTFQFGQFVRPVCLPEPGEHFNAGFICTAGWGRLSEGGRLPQVLQQVNLPILTQEECEAV
 LLTLKNPITGKTFCTGSPDGGRDACQDSSGSLMCQNRKGAWTLAGVTSWGLGGRSWRNARKKEQGS
 PGIFTDLRRVLPWILKHIQTGHRKSTKALCSEPDGLISGSEGLHFPESLHLYYESKQLCVWTFVLPED
 MHMLNLNLSHLDIESCHHNYLAMYSLEDRLVGKFCGENLPSSILIGSSSIRLRFISDATDYATGFNLTYKA
 LKPRYHPDSCGRSLTILFEEGTIQSLHYPEDYSNMASCTWIFQAPNYCIKLSFQSLIEENGDCSSDYV
 TVHSDVEKEKEIARFCDYVIPSPVLSSSSVMLISFQSDENGARGFQADISFISRADLNISISEDESVP
 ETWDLPPGAMEIFDAERDTHTKPPYEEDIGEMPAIDSGLLKQGERRGKH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9082_f08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_172908

ORF Size: 1827 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_172908.3](#), [NP_766496.2](#)

RefSeq Size: 1924 bp

RefSeq ORF: 1830 bp
 Locus ID: 244199
 UniProt ID: [Q7M761](#)
 Cytogenetics: 7 E3
 MW: 67.6 kDa

Product images:



Circular map for MR223336