

## Product datasheet for **MR223379L4V**

### Dok5 (NM\_001163686) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dok5 (NM_001163686) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dok5
Synonyms:	2700055C10Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001163686
ORF Size:	912 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223379).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001163686.1</a> , <a href="#">NP_001157158.1</a>
RefSeq Size:	1768 bp
RefSeq ORF:	915 bp
Locus ID:	76829
UniProt ID:	<a href="#">Q91ZM9</a>
Cytogenetics:	2 92.26 cM



[View online »](#)

**Gene Summary:**

DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK5 functions in RET-mediated neurite outgrowth and plays a positive role in activation of the MAP kinase pathway. Putative link with downstream effectors of RET in neuronal differentiation. [UniProtKB/Swiss-Prot Function]