

Product datasheet for **MR223383**

Ablim2 (NM_177678) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ablim2 (NM_177678) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ablim2
Synonyms:	A1606905; C230091L11
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR223383 representing NM_177678
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCGCAGTGTGCGAGCCCCAGGCTGCTCATGCCCACTGGAGAAGCCGGCCAGCACTGCGATCCTGT
 GTAATACCTGTGGGAATGTGTGCAAGGGAGAGGTGCTGCGTGTGCAAGAACAAGTACTTCCACATCAGGTG
 CTTTGTGTGCAAAGCGTGTGGCTGTGACCTGGCTGAGGGTGGCTTCTTCGTGAGGCAGGGTGAGCACATT
 TGCACACGTGACTACCAGCGGCTCTATGGCACTCGCTGCTTCAGCTGTGACCGCTTCATTGAAGCGAAG
 TGGTGTGCGGCGCTTGGGAAGACCTACCACCCTGACTGCTTCGTGTGTGCTGTCTGCAAGGCTGCCCTTTC
 TCCTGGGACCGTGTGACCTTCAACGGGAAGGAGTGTATGTGCCAGAAGTGTCCCAACCCACATTGCTG
 GGCAACAGTGTCTACGTGGCCAGGGCTCCGGAGTTGTGGGGCTGTGGCTTAGAAATCAAGAATGGCC
 AGGCCCTGGTGGCTTTGGACAAGCACTGGCACCTGGGATGCTTCAAGTGAAGACTTGTGGGAAACTCCT
 CAATGCAGAGTACATCAGCAAGGACGGGCTGCCCTACTGCGAGGCTGATTATCACAGCAAGTTTGGCATC
 CGCTGCGACGGCTGTGAGAAGTACATCACGGGCCGGGTGCTAGAGGCCGGAGAGAAGCACTACCACCTT
 CATGTGGCTATGCGTCAGGTGCGGCCAGATGTTCTCAGAGGGCGAGGAGATGTACCTGCAAGGCTCCTC
 CATCTGGCATCCGGCATGTGACAGGCAGCCAGGACTGAAGACAAGAGCAAGGAAACCAGGACTTCTTCG
 GAGAGCATTGTCTCCGTGCCTCCAGCACGTACGGTCCCGGAGCCGTGTGATCTACGCAAAGCTTG
 GCGATGAGATCCTGGACTACAGGACTTGGCTGCTTCCCAAAAACAAGCGATCTACAACATCGACCG
 CCCCAGATGATCCTATTGCCCCTACATTAGCCACTCGGCCGTGGGGACAGGCAGAGCTATGGCGAG
 GGGGATCAGGATGATCGGTCTACAAGCAGTGCAGGACCTCCAGCCCCAGCTCTGCCGGCTCAGTCAGCC
 TCGGACACTACCCCCAACCTCACGGTCACCCCAGCACTACAGTGTCCAGGCAGCGAAAGTGGCCGTAG
 CACCCCGAGTCTCTCAGTGCCTCCGACAGCAGACCTCCTTCTCTACCTACCAGCAGGCTCCGAGACAC
 TTCCACGTCCCAGACACTGGCGTAAAAGATAACATCTATAGGAAACCCCTATCTACAAACAGCAGGCTG
 CCAGGCGATTGGACGTGGAGGACAGCAGCTTTGACCAGGACAGCAGGAAGAAGACAACCTGGCTGCTTCT
 CAAGGGGGATGCAGACACCAGAATAACTCTCCAGACCTGGACAGTCAAGTCGCTGTCCCTCAGCAGCGGG
 ACTGACCAGGAACCTCTCAAAGAATGGCAGGGGACAGCCTCTACTCACGCTTTCCTTATCCAAACCTG
 AACTCTTCCAGGACCCAGGAAGGATGGCTGGACCTCCGGAATGCCAACCTGGCCCTTGTGGAGCAGA
 CCCGGATGCCAGTGGGGCACGCGAGAGTATAAGATATATCCCTATGATTCCTCATCGTAACAAACCGA
 ATTCGTGTGAAACTGCCAAAGACGTGGACCGGACACGGCTGGAGAGACCTATCACCTGAGGAGTTCC
 AAGAAAGTGTGGGATGAGCATCGAGGAGTTTGACCGCCTGGCCCTGTGGAAGAGGAATGACCTCAAGAA
 GAAAGCCTTGCTGTTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR223383 representing NM_177678
 Red=Cloning site Green=Tags(s)

MSAVSQPQAHAHAPLEKPASTAILCNTCGNVCKGEVLRVQNKYFHIRCFVCKACGCDLAEGGFFVRQGEHI
 CTRDYQRLYGTRCFSCDRFIEGEVVSALGKTYHPDCFVCAVCRLPFPDPGDRVTFNGKECMQKCSPTLL
 GNSAHVAQGLRSCGGCGLEIKNGQALVALDKHWHLGCFKCKTCGKLLNAEYISKDGLPYCEADYHSKFGI
 RCDGCEKYITGRVLEAGEKHYPSCALCVRGQMFSEGEEMYLQGSIIWHPACRQAARTEDKSKETRSTSS
 ESIVSVPASSTSGSPSRVIYAKLGDEILDYRDLAALPKNKAIYNIDRPMISYSPYISHSAVGDRQSYGE
 GDQDDRSYKQCRSTSSPSSAGSVSLGHYTPTSRSPQHYSRPGSESGRSTPSLSVHSDSRPPSSTYQQAPRH
 FHVPDTPGKDNIRKPPYKQHAARRLDVEDSSFDQDSRKKTTWLLKGDADTRTNSPDLDSQSLSSG
 TDQEPLQRMAGDSLYSRFPYKPDTPGPRKDGDLRNLNANLAPCGADPDASWGTREYKIYPYDSLIVTNR
 IRVKLPKDVDRTRLERHLSPEEFQEVFGMSIEEFDRLALWKRNDLKKKALLF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_177678

ORF Size: 1836 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_177678.7 , NP_808346.3
RefSeq Size:	3570 bp
RefSeq ORF:	1839 bp
Locus ID:	231148
UniProt ID:	Q8BL65
Cytogenetics:	5 B3
MW:	68.6 kDa
Gene Summary:	May act as scaffold protein. May stimulate ABRA activity and ABRA-dependent SRF transcriptional activity.[UniProtKB/Swiss-Prot Function]