

Product datasheet for **MR223397L3V**

Arhgef7 (NM_017402) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgef7 (NM_017402) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgef7
Synonyms:	betaPix; betaPix-b; betaPix-c; Cool; cool-1; mKIAA0142; p85Cool1; p85SPR; Pak3bp; PIX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017402
ORF Size:	1938 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017402.4
RefSeq Size:	4488 bp
RefSeq ORF:	1941 bp
Locus ID:	54126
UniProt ID:	Q9ES28
Cytogenetics:	8 A1.1



[View online »](#)

Gene Summary:

Acts as a RAC1 guanine nucleotide exchange factor (GEF) and can induce membrane ruffling. May function as a positive regulator of apoptosis. Functions in cell migration, attachment and cell spreading. Promotes targeting of RAC1 to focal adhesions. Downstream of NMDA receptors and CaMKK-CaMK1 signaling cascade, promotes the formation of spines and synapses in hippocampal neurons (By similarity).[UniProtKB/Swiss-Prot Function]