

## Product datasheet for **MR223406L4V**

### Gtf2a1 (NM\_031391) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gtf2a1 (NM_031391) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gtf2a1
Synonyms:	6330549H03Rik; AA536742; AA959775; AW060250; Tfii1; TfiiAa/b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031391
ORF Size:	1134 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223406).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031391.2</a> , <a href="#">NP_113568.2</a>
RefSeq Size:	5704 bp
RefSeq ORF:	1137 bp
Locus ID:	83602
UniProt ID:	<a href="#">Q99PM3</a>
Cytogenetics:	12
Gene Summary:	TFIIA is a component of the transcription machinery of RNA polymerase II and plays an important role in transcriptional activation. TFIIA in a complex with TBP mediates transcriptional activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)