

Protein Sequence:	>MR223412 representing NM_024469 Red=Cloning site Green=Tags(s) MDEGIPHLQERQLLEHRDFIGLDYSSL YMCKPKRSLKRDDTKDITYKLPHRLIEKKRRDRINECIAQLKDL LPEHLKLTTLGHLEKAVVLEL TLKHLKALTAL TEQQHQKIIALQNGERSLKSPVQADLDAFHSGFQTCAK EVLQYLARFESWTPREPRCAQLVSHLHAVATQLLTPQVPSGRGSGRAPCSAGAAAASGPervARCVPIQ RTQPGTEPEHDTDTDSGYGGEAEQGRAAVKQEP PGDSSAPKRPKLEARGALLGPEPALLGSLVALGGA PFAQPAAPFCLPFYLLSPSAAAYVQPWLDKSGLDKYL YPAAAPFPLLYPGIPAAAAAAAAAFAFCLSS VLSPPPEKAGATAGAPFLAHEVAPPGLRPQHAHSRTHLPRAVNPESQEDATQPAKDAP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_024469
ORF Size:	1230 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024469.1</u> , <u>NP_077789.1</u>
RefSeq Size:	1421 bp
RefSeq ORF:	1233 bp
Locus ID:	79362
UniProt ID:	<u>Q99PV5</u>
Cytogenetics:	6 77.7 cM
MW:	44.4 kDa