

Product datasheet for **MR223429L3V**

Glmn (NM_001161738) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Glmn (NM_001161738) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Glmn
Synonyms:	9330160J16Rik; AW227515; Fap48; Fap68
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001161738
ORF Size:	1788 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223429).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001161738.1 , NP_001155210.1
RefSeq Size:	2001 bp
RefSeq ORF:	1791 bp
Locus ID:	170823
UniProt ID:	Q8BZM1
Cytogenetics:	5 E5



[View online »](#)

Gene Summary:

Regulatory component of cullin-RING-based SCF (SKP1-Cullin-F-box protein) E3 ubiquitin-protein ligase complexes. Inhibits E3 ubiquitin ligase activity by binding to the RING domain of RBX1 and inhibiting its interaction with the E2 ubiquitin-conjugating enzyme CDC34. Inhibits RBX1-mediated neddylation of CUL1 (By similarity). Required for normal stability and normal cellular levels of key components of SCF ubiquitin ligase complexes, including FBXW7, RBX1, CUL1, CUL2, CUL3, CUL4A, and thereby contributes to the regulation of CCNE1 and MYC levels (PubMed:22405651). Essential for normal development of the vasculature (PubMed:22405651). Contributes to the regulation of RPS6KB1 phosphorylation (By similarity).[UniProtKB/Swiss-Prot Function]