

## Product datasheet for **MR223445L3V**

### Adamdec1 (NM\_021475) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adamdec1 (NM_021475) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adamdec1
Synonyms:	2210414L24Rik; AI574231; Dcsn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021475
ORF Size:	1401 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223445).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021475.2</a> , <a href="#">NP_067450.1</a>
RefSeq Size:	2412 bp
RefSeq ORF:	1404 bp
Locus ID:	58860
UniProt ID:	<a href="#">Q9R0X2</a>
Cytogenetics:	14 D1
Gene Summary:	May play an important role in the control of the immune response and during pregnancy. [UniProtKB/Swiss-Prot Function]



[View online »](#)