

## Product datasheet for **MR223452L1V**

### Tnfsf4 (NM\_009452) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnfsf4 (NM_009452) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnfsf4
Synonyms:	Ath-1; Ath1; CD134L; gp34; OX-40L; Ox40l; Tnlg2b; TXGP1; Txgp1l
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_009452
ORF Size:	594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223452).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009452.2</a> , <a href="#">NP_033478.1</a>
RefSeq Size:	1609 bp
RefSeq ORF:	597 bp
Locus ID:	22164
UniProt ID:	<a href="#">P43488</a>
Cytogenetics:	1 69.75 cM
Gene Summary:	Cytokine that binds to TNFRSF4. Co-stimulates T-cell proliferation and cytokine production. [UniProtKB/Swiss-Prot Function]



[View online »](#)