

Product datasheet for **MR223478**

Hmga1 (NM_001166543) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hmga1 (NM_001166543) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Hmga1
Synonyms: AL023995; Hmg; HMG-I(Y; Hmga; Hmga1a; Hmga1b; Hmgj; HMGI(Y; Hmgij; Hmgy
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR223478 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCGAGTCGGGCTCAAAGTCCAGCCAGCCCTGGCCTCAAGCAGGAAAAGGATGGGACTGAGAAGC
GAGGCCGGGCGAGGCCACGCAAGCAGCCTCCGGTGAGTCCTGGGACGGCGCTGGTCGGGAGTCAGAAAGA
GCCAGTGAAGTGCCAACTCCGAAGAGACCTCGGGCCGACCAAAGGGAAGCAAGAATAAGGGCGCCGCC
AAGACCCGAAAAGTCACCACAGCTCCAGGGAGGAAACCAAGGGGCAGACCAAGAAACTGAAGGAGGAAG
AGGAGGGCATCTCCAGGAGTCCTCTGAGGAGGAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223478 protein sequence
Red=Cloning site Green=Tags(s)

MSESGSKSSQPLASKQEKDGTETKRGRRPRKQPPVSPGTALVGSQKEPSEVPTPKRPRGRPKGSKNKGAA
KTRKVTTAPGRKPRGRPKKLKEEEEGISQESSEEEQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001166543

ORF Size:

318 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166543.1</u> , <u>NP_001160015.1</u>
RefSeq Size:	1714 bp
RefSeq ORF:	321 bp
Locus ID:	15361
Cytogenetics:	17 14.5 cM
MW:	11.5 kDa
Gene Summary:	This locus encodes a member of the nuclear, non-histone high mobility group protein family. This architectural transcription factor binds to A-T rich DNA sequences and participates in enhanceosome formation, chromatin remodeling and regulation of transcription. This protein functions in many cellular processes, including cell growth and differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Oct 2009]