

Product datasheet for **MR223488A1V**

Mouse Tor1aip2 (NM_001160180) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Tor1aip2 (NM_001160180) AAV Particle
Tag: Myc-DDK
Symbol: Tor1aip2
Synonyms: 15kDa; 1110020D10Rik; A130072J07; AA103493; AW060462; AW610675; C77739; Ifrg15; Lull1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR223488 representing NM_001160180
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTCTCCGATAATTCACACTGCCCTGATTGTGGCAGCAGTGGTTCCTAGTTTGAAGTACTAGGGCACT
 GGTTGTACAAACTGAACCTGTTGAAAATGAATGTTACCAAGTGTCTTAGACCGTATTAACAGGGCTGA
 CTACTGCCCTGAGTGTACCCCTGACAATCCTGCTAATAGAAGCCTTGTCTGCCTGGTCTTTCCCACTT
 GAGTGGGCACCTCAAAATCTTACCAGGTGGACCTTTGAAAAGCTTGCCATCCATTTCTCTGGGTCCTC
 CACTGGTTAGAAAAAGATACATGACTCCAGGGTAGCTGGCTTAACTCCTGCATTACAATTAATTTGTG
 CAGAACAGACAAAACCTTAAACAAGAACTTGGCCAAAGCAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223488 representing NM_001160180
 Red=Cloning site Green=Tags(s)

MFSDNSHCPDCGQWFPSLELGHWLYQTELVENECYQVFLDRINRADYCPICYPDNPANRSLVLPWSFPL
 EWAPQNLTRWTFEKACHPFLGPPVLRKKIHDSRVAGFNPALQLILSRDKTLNKKLGQSK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001160180



[View online »](#)

ORF Size:	396 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001160180.1</u> , <u>NP_001153652.1</u>
RefSeq Size:	2697 bp
RefSeq ORF:	396 bp
Locus ID:	240832
UniProt ID:	<u>Q9ER81</u>
Cytogenetics:	1 G3
MW:	15.3 kDa