

## Product datasheet for **MR223488L3V**

### Tor1aip2 (NM\_001160180) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tor1aip2 (NM_001160180) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tor1aip2
Synonyms:	15kDa; 1110020D10Rik; A130072J07; AA103493; AW060462; AW610675; C77739; Ifrg15; Lull1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001160180
ORF Size:	396 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223488).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001160180.1</a> , <a href="#">NP_001153652.1</a>
RefSeq Size:	2697 bp
RefSeq ORF:	396 bp
Locus ID:	240832
UniProt ID:	<a href="#">Q9ER81</a>
Cytogenetics:	1 G3
Gene Summary:	Required for endoplasmic reticulum integrity. Regulates the distribution of TOR1A between the endoplasmic reticulum and the nuclear envelope as well as induces TOR1A, TOR1B and TOR3A ATPase activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)