

## Product datasheet for **MR223504L3V**

### Rnf152 (NM\_178779) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnf152 (NM_178779) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rnf152
Synonyms:	A930029B02Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178779
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223504).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178779.4</a> , <a href="#">NP_848894.1</a>
RefSeq Size:	8529 bp
RefSeq ORF:	612 bp
Locus ID:	320311
UniProt ID:	<a href="#">Q8BG47</a>
Cytogenetics:	1 E2.1


[View online »](#)

**Gene Summary:**

E3 ubiquitin-protein ligase mediating 'Lys-63'-linked polyubiquitination of RRAGA in response to amino acid starvation. Thereby, regulates mTORC1 signaling and plays a role in the cellular response to amino acid availability (PubMed:25936802). Also mediates 'Lys-48'-linked polyubiquitination of target proteins and their subsequent targeting to the proteasome for degradation. Induces apoptosis when overexpressed (By similarity).[UniProtKB/Swiss-Prot Function]