

Product datasheet for **MR223537L4V**

Pgap1 (NM_001163314) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pgap1 (NM_001163314) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pgap1
Synonyms:	5033403E17Rik; 9030223K07Rik; A530084K22; D230012E17Rik; oto
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001163314
ORF Size:	2766 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223537).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001163314.2 , NP_001156786.1
RefSeq Size:	2871 bp
RefSeq ORF:	2769 bp
Locus ID:	241062
UniProt ID:	Q3UUQ7
Cytogenetics:	1 C1.2
Gene Summary:	Involved in inositol deacylation of GPI-anchored proteins. GPI inositol deacylation may important for efficient transport of GPI-anchored proteins from the endoplasmic reticulum to the Golgi (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)