

## Product datasheet for **MR223547L3V**

### **Nrp2 (NM\_001077407) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Nrp2 (NM_001077407) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nrp2
Synonyms:	1110048P06Rik; Np-2; Np2; Npn-2; Npn2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001077407
ORF Size:	2703 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223547).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001077407.1</a> , <a href="#">NP_001070875.1</a>
RefSeq Size:	4696 bp
RefSeq ORF:	2706 bp
Locus ID:	18187
Cytogenetics:	1 C2
Gene Summary:	High affinity receptor for semaphorins 3C, 3F, VEGF-165 and VEGF-145 isoforms of VEGF, and the PLGF-2 isoform of PGF.[UniProtKB/Swiss-Prot Function]



[View online »](#)