

Product datasheet for **MR223604L4V**

Nnat (NM_180960) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nnat (NM_180960) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nnat
Synonyms:	5730414I02Rik; AW107673; Peg5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_180960
ORF Size:	162 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223604).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_180960.2
RefSeq Size:	1142 bp
RefSeq ORF:	165 bp
Locus ID:	18111
UniProt ID:	Q61979
Cytogenetics:	2 H1
Gene Summary:	May participate in the maintenance of segment identity in the hindbrain and pituitary development, and maturation or maintenance of the overall structure of the nervous system. [UniProtKB/Swiss-Prot Function]



[View online »](#)