

## Product datasheet for **MR223606**

### Ube3a (NM\_173010) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ube3a (NM_173010) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ube3a
Synonyms:	4732496B02; 5830462N02Rik; A130086L21Rik; Hpve6a; mKIAA4216
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide  
Sequence:**

>MR223606 representing NM\_173010  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAAGCGAGCAGCTGCAAAGCATCTAATAGAAGCTACTACCATCAGTTAACTGAGGGCTGTGGAATG  
 AGGCCTGCACGAATGAGTTTTGTCTTCTGTCCAACCTTTCTTCGTATGGATAACAATGCAGCAGCTAT  
 TAAAGCCCTTGAGCTTTATAAAATTAATGCAAACTCTGTGATCCTCATCCCTCCAAGAAAGGAGCAAGC  
 TCAGCTTACCTTGAGAACTCAAAGGTGCATCTAACTCAGAGATAAAAATGAACAAGAAGGAAGGAA  
 AAGATTTAAAGATGTGATTTACCTAACTGAAGAGAAAGTATATGAAATTTATGAATTTGTAGAGAGAG  
 TGAGGATTATCCCTTTAATTCGTGTAATTGGAAGAATATTTCTAGTGCTGAGGCACTGGTTCTGAGC  
 TTTCCGAAAGTCAAACAGCACACAAAGGAGGAATTGAAATCTCTCAAGAAAAGGATGAAGACAAGGATG  
 AAGATGAAAAGGAAAAAGCTGCATGTTCTGCTGCTGCTATGGAAGAAGACTCAGAAGCATCTTCTCAAG  
 GATGGGTGATAGTTCACAGGGAGACAACAATGTACAAAATTAGGTCCTGATGATGTGACTGTGGATATT  
 GATGCTATTAGAAGGGTCTACAGCAGTTTGTCTCGCTAATGAAAAATTAGAACTGCCTTCTGAATGCAC  
 TTGTATATCTGTACCTAACGTGGAATGTGATTTGACATATCATAATGTGTATACTCGAGATCCTAATTA  
 TCTCAATTTGTTCAATTATTGTAATGGAGAATAGTAATCTCCACAGTCTGAAATATCTGGAATGGCGTTG  
 CCATTATTTGCAAAGCTATGTGTAAGCTACCCCTTGAAGCTCAAGGAAAAGTATTAGGCTGTGGTCTA  
 AATACAGTGCTGACCAGATTCGGAGAATGATGGAACATTTTCAAGCACTTATTACCTACAAAAGTCATAAG  
 CAATGAATTTAATAGCCGAAATCTAGTGAATGATGATGATGCCATTGTTGCTGCTCAAAGTGTGTTGAAA  
 ATGGTTTACTATGCAAAATGATGAGGGAGGGGATGTGGACAAAATCATAATGAGGAAGATGATGAAGA  
 CCATACCTGAGTCCAGCGAATTAACACTTCAGGAGCTTCTGGGAGATGAAAGAAGAAATAAGAAAGGTCC  
 TCGAGTGGATCCACTAGAAACCGAACTTGGCGTTAAAACCTAGACTGTGAAAACCACTTATCTCCTTT  
 GAAGAATTCATTAAATGAACCACTGAATGATGTTCTAGAAATGGACAAAAGATTATACCTTTTTCAAAGTTG  
 AAACAGAGAACAATTTCTTTTTATGACATGTCCCTTTATATTGAATGCTGTACAAAAGAATCTGGGATT  
 ATATTATGACAATAGAATTCGCATGTACAGTGAAGAAGAATCACTGTTCTTTACAGCCTAGTTCAAGGA  
 CAGCAGTTGAATCCGATTTTGGACTCAAAGTCAAGCTGACCATATTATAGATGATGCACTGGTCCGGC  
 TAGAGATGATTGCTATGGAAAATCCTGCAGACTGAAGAAGCAGTTGTATGTGGAATTTGAAGGAGAACA  
 AGGAGTAGATGAGGGAGGCGTTTCAAAGAGTTTTTTTTCAGTTGGTTGTGGAGGAAATTTTAAATCCAGAT  
 ATTTGGTATGTTACATATGATGAAGCTACGAAATTTTGGTTTAAATCCATCTTCTTTGAAACTGAGG  
 GTCAGTTTACTCTGATTGGCATAGTCTGGGTCTGGCTATTTACAATAATTGTACTGGATGTCCATTT  
 TCCCATGGTTGTATACAGGAAGCTAATGGGGAAAAAGGAACCTTTCTGTGACTTGGGAGACTCTACCCCA  
 GTTTTATATCAGAGTTTAAAGGATTTATTGGAATATGAAGGGAGTGTGGAAGATGATATGATGATCACTT  
 TCCAGATATCACAGACAGATCTTTTTGGTAACCAATGATGTATGATCTAAAAGAAAATGGTGATAAAAT  
 TCCAATTACAAATGAAAACAGGAAGGAATTTGTCAATCTCTATTCAAGACTACATTCTCAATAAATCTGTA  
 GAAAAACAATTCAGGCATTTTCGAGAGGTTTTTCATATGGTGACTAATGAATCGCCCTTAAAATACTTAT  
 TCAGACCAGAAGAAATGAATTGCTTATATGTGGAAGCCGGAATCTAGATTTCCAGGCACTAGAAGAAAC  
 TACAGAGTATGACGGTGGCTATACGAGGGAATCTGTTGTGATTAGG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


**Protein Sequence:** >MR223606 representing NM\_173010  
Red=Cloning site Green=Tags(s)

MKRAAAKHLIERYYHQLTEGCGNEACTNEFCASCPTFLRMDNNAAAIKALELYKINAKLCDPHPSKKGAS  
 SAYLENSKGASNNSEIKMNKKEGKDFKDVLYL TEEKVYEIYEFCRESEDYSPLIRVIGRIFSSAEALVLS  
 FRKVKQHTKEELKSLQEKDEDEDEDEKEKAACSAAMEEDSEASSRMGDSSQGDNNVQKLGPDVTVDI  
 DAIRRVYSLLANEKLETAFLNALVYLSPNVECDLTYHNVTYTRDPNYLNLFIIVMENSNLHSPEYLEMAL  
 PLFCAMCKLPLEAQGKLIRLWSKYSADQIRRMETFQQLITYKVISNEFNSRNLVDDDAIVAASKCLK  
 MVYYANVVGDDVDTNHNEEDDEEPIPESEL TLQELLGDERRNKKGPRVDPLETELGVKTLDCRKPLISF  
 EEFINEPLNDVLEMDKDYTFK VETENKFSFMTCPFILNAVTKNLGLYYDNRIRMY SERRITVLVLSLVQG  
 QQLNPYLRLKVRDHIIDDALVRLEMIAMENPADLKKQLYVEFEGEQGVDEGGVSKFEFFLVVEEIFNPD  
 IGMFTYDEATKLFWFNPSSFETEGQFTLIGIVLGLAIYNNCILDVHFPVMVYRKL MGKKGTRFDLGDSPH  
 VLYQSLKDLLEYEGSVEDDMITFQISQTDLFGNPMMYDLKENGDKIPITNENRKEFVNL YSDYILNKSV  
 EKQFKAFRRRGFHMVTNESPLKYLFRPEEIELLICGSRNLDFQALEETTEYDGGYTRESVVIR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_173010

**ORF Size:** 2286 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_173010.3</a> , <a href="#">NP_766598.1</a>
<b>RefSeq Size:</b>	3888 bp
<b>RefSeq ORF:</b>	2289 bp
<b>Locus ID:</b>	22215
<b>UniProt ID:</b>	<a href="#">O08759</a>
<b>Cytogenetics:</b>	7 33.95 cM
<b>MW:</b>	88.1 kDa
<b>Gene Summary:</b>	<p>E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and transfers it to its substrates. Several substrates have been identified including the ARNTL/BMAL1, ARC, RAD23A and RAD23B, MCM7 (which is involved in DNA replication), annexin A1, the PML tumor suppressor, and the cell cycle regulator CDKN1B (PubMed:20211139, PubMed:24728990). Additionally, may function as a cellular quality control ubiquitin ligase by helping the degradation of the cytoplasmic misfolded proteins. Finally, UBE3A also promotes its own degradation in vivo (By similarity). Plays an important role in the regulation of the circadian clock: involved in the ubiquitination of the core clock component ARNTL/BMAL1, leading to its proteasomal degradation (PubMed:24728990). Acts as a regulator of synaptic development by mediating ubiquitination and degradation of ARC (PubMed:20211139). Synergizes with WBP2 in enhancing PGR activity (By similarity). [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR223606