

Product datasheet for **MR223607L3V**

Naip6 (NM_010871) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Naip6 (NM_010871) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Naip6
Synonyms:	Birc1f; Naip-rs4; Naip-rs4A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010871
ORF Size:	4209 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223607).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010871.2 , NP_035001.2
RefSeq Size:	6705 bp
RefSeq ORF:	4212 bp
Locus ID:	17952
UniProt ID:	Q9JIB6
Cytogenetics:	13 53.09 cM
Gene Summary:	Closest sequence match is AF381772. [provided by RefSeq, Jul 2008]



[View online »](#)