

Product datasheet for **MR223627**

H3f3a (NM_008210) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H3f3a (NM_008210) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: H3f3a
Synonyms: EyeLinc14; H3-3; H3-3a; H3-3b; H3.3; H3.3A
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR223627 representing NM_008210
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCGTACAAAGCAGACTGCCCGCAAATCCACCGGTGGTAAAGCACCCAGGAAACAACCTGGCTACAA
AAGCCGCTCGCAAGAGTGCGCCCTCTACTGGAGGGGTGAAGAAACCTCATCGTTACAGGCCTGGTACTGT
GGCCCTCCGTGAAATCAGACGCTATCAGAAGTCCACTGAACCTTCTGATCCGCAAGCTCCCTTTCCAGCGT
CTGGTGCAGAAATTGCTCAGGACTTCAAACAGATCTGCGCTTCCAGAGTGCAGCTATTGGTGCTTTGC
AGGAGGCAAGTGAGGCCTATCTGGTTGGCCTTTTTGAAGATACCAATCTGTGTGCTATCCATGCCAAACG
TGTAACAATTATGCCAAAAGATATCCAGCTAGCACGCCGCATACGCGGAGAACGTGCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223627 representing NM_008210
Red=Cloning site Green=Tags(s)

MARTKQTARKSTGGKAPRKQLATKAARKSAPSTGGVKKPHRYRPGTVALREIRRYQKSTELLIRKLPFQR
LVREIAQDFKTLRFQSAALQEAAYLVGLFEDTNLCAIHAKRVTIMPKDIQLARRIRGERA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9034_d03.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_008210

ORF Size: 408 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008210.5, NP_032236.1](#)
RefSeq Size: 1076 bp
RefSeq ORF: 411 bp
Locus ID: 15078
UniProt ID: [P84244](#)
Cytogenetics: 1 H4
MW: 15.8 kDa

Gene Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene contains introns and its mRNA is polyadenylated, unlike most histone genes. The protein encoded is a replication-independent member of the histone H3 family. [provided by RefSeq, Nov 2015]

Product images:



Circular map for MR223627