

Product datasheet for **MR223648**

G2e3 (NM_001167964) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	G2e3 (NM_001167964) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	G2e3
Synonyms:	6030408C04Rik; 9030416F18; AW046379; D930034K21Rik; mKIAA1333
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR223648 representing NM_001167964
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAATGAAAATAAACCTGATAATTCACAAAGTCTTGCTTGTGTTTTCTGCCAAAAAATGATGACTGTC
 CCAATAAATACGAGAGAAAAGAAAACATATGAAAAATGGAATTTTGTGTCCTACTACTGCTTGTAAAT
 GTCAAGTGGAAATTTGGCAGAGAGGTAAGAAGAAGAAGGAGTTTATGGCTTCTAATAGAAGACATCAGG
 AAGGAAGTCCAGAGAGCCTCAAACGAAGTGTACAGTTTGAAGAAGAATGGTCTTCAATTGGATGTG
 TGGTACCCACGTGTAACGAAGTTATCACTTGCCTTGTGGACTTCAGAAAGAATGATTTTCCAGTTTAC
 TGACAATTTTGCCTCATTTTGTGAAACATCGACCCGTTCAAGCAATTACATCTAATAAATATAGTAGT
 TCTTTACCATGTACTATCTGTTTGGAAATTTGTCGAGCCTATTCTACTTATAACATACTACAAGTCCCT
 GTTGTAAAGATGCTTGGTTTCATAGAGACTGTTTACAGGTCCAAGCAATAAATGCAGGAGTGTTCCTT
 TAGATGTACACTATGCAATAATACAGACATCTTTCAGAAGGAGATGTTGAGAATGGGAATTCATATTCCT
 GAAAAAGATGCATCCTGGGAATTAGAAGAAAACGCGTATCAAGAGCTGCTGCAGAGCCACGACCGCTGTG
 ATATTCGAAGATGTCCTGCAAAAAAGGCCGAGACTACAACGAACCTAACAGCAATGGGAAGTAAAGCG
 GTGTCAGAGTTGTGGTTCTAGTGAACACACTTACGCTGTTTCTCACTACAGTCTTGGGAACAAAATTGG
 GAGTGTGGACTGTAGAAGAATTACCTACACTTACAGATTTCCAAAAGGCCCCCAACACCCATTAGCCA
 ACTCTACTAATGTGACAGTTACTGACTGCTTGTGGAAGATCATCATCTAAGCTCCCGAGACAGTCAAC
 TGTAGCCAGCATAAAGAAGTACTGAGTTATTTTATCTGTATTTACAGCTGTGAAAGAGAATCTTTAC
 TATGAAGTTGGCAAAATGCTTGCCATTTCTTTAGTTTATGTTGTTGTTGCTTTCGCTGGTTTCTTTTCTGAAA
 CCTTGTAACTGCCTTGTATGACCAGAAAAACTCTACCAACTTTGGATGATGATCAGACATTGA
 TGTGACAGATTATAATCAAGATAGATTCTGCAACTGATTTAAACATTTTAAATTCAGTAATAAGCCAA
 CACTATAACTACTTAGAGGTCAGTGGATGTTTAAAGACTTACGACATCATTAAAGTATAAATTCATGTTGG
 TAAAAGACATCCTTTTCTACCATGTAATTAATAGAGTCAAGGCACCTTTTGGAGTTTTAAACAGGGTCT
 GAAAACGCTTGGTGTGTTTAGAGAAAATCAAACCTACCCAGAAGCATTATAAGATCCTATGTCATAAA
 CCTGAGAATCTCTGCAAAAAACCTCAGTACCTTTCACAATTCACCTCGGTAGCGGACGTACAGACTT
 TGAGGTTTTGGAACAGTTACTTAAAGGCTATAGAGGATGGTAAATCTGCAACAACGATGGAAGATATTCT
 TATTTTGAACCTGGCTGTAGTCTGTCCCAACTGGATTTAAACCCTCTCTGTCAGTTGAATGTCTA
 CATGTGGATTTTCTGTTGCAGACAAATATAGAAACATTTAGTTCTACCAGCCACCAACACGTATGAAG
 AGTTCCAGAAAATATGGATTTACCATAAGAGACACTTAAGATTAGAAAAGGAAGAAAGGTCTCACAT
 ACTACCTCGGACATTAATGTTTCTTCTAATGAAGAAATGCTTATT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR223648 representing NM_001167964
 Red=Cloning site Green=Tags(s)

MNENKPDNSQSLACVFCRKNDDCPNKYGEKTYEKWVSVHYCLLMSSGIWQRGKEEYVYGFLLIEDIR
 KEVQRASKLKCTVCKKNGASIGVVPCKRSYHLPCGLQKECIFQFTDNFASFCKWHRPVQAITSNKYSS
 SLPCTICLEFVEPIPTYNILQSPCKNAWFHRDCLQVQAINAGVFFRCTLNNTDIFQKELRMGIHIP
 EKDASWELEENAYQELLQSHDRCDIRRHCKKGRDYNENPSKWEVKRCQSCGSSGTHLACSSLQSWEQNW
 ECLDCRRITYTSDFQKAPKHLANSTNVTVDCLLEESSKLRQSTVAQHKELLSLFLSVFTAVKENLY
 YEVGKMLAISLVHGGPSPGFFSETLFNCLAYGPENTLPTLDDVSDIDVAQIIKIDSATDLNINLSVISQ
 HYNYLEVSGCLRLLTSLSDKFMLVKDILFYHVINRVKAPFESFKQGLKTLGVLEKIQTTYPEAFYKILCHK
 PENLSAKNLSLFTIHSVADVQTLRFWNSYLKAIEDGKSATTMEDILIFATGCSSVPPTGFKPSLSVECL
 HVDFPVADKYRNHLVLPATNTYEEFQENMDFIIRDTLRLEKEERSHILPRTLNVSSNEEMLI


SGP**TRTRRLEQKLI**SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001167964

ORF Size: 1866 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001167964.1 , NP_001161436.1
RefSeq Size:	6352 bp
RefSeq ORF:	1869 bp
Locus ID:	217558
UniProt ID:	Q5RIY2
Cytogenetics:	12 B3
MW:	71.5 kDa
Gene Summary:	E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates (By similarity). Required for prevention of apoptotic death in early embryogenesis. [UniProtKB/Swiss-Prot Function]