

Product datasheet for **MR223656**

Satb1 (NM_001163631) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Satb1 (NM_001163631) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Satb1
Synonyms:	2610306G12Rik; AW413156
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR223656 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATCATTTGAACGAGGCAACTCAGGGGAAAGAACATTCAGAAATGTCTAACAATGTGAGTGATCCGA
 AGGGTCCACCCGCAAGATTGCCCGCTGGAGCAGAACGGGAGCCCTCTAGGAAGAGGAAGGCTTGGGAG
 CACAGGTGAAAAGATGCAGGGAGTGCCTTTAAACACTCGGGCCATCTCATGAAAACCAACCTTAGGAAA
 GGAACCATGTTACCAGTTTTCTGCGTGGTGGAACATTATGAAAACGCCATTGAGTATGATTGCAAGGAGG
 AGCACGCGGAATTTGATTGGTGAGAAAGGATATGCTTTTCAACCAGCTGATAGAGATGGCGTTGCTGTC
 TCTAGGCTATTCACACAGCTCTGCTGCCAAGCAAAGGGCTCATCCAGGTTGGGAAGTGAATCCAGTT
 CCACTGTCGATGTGACAGATGCCCTGATGCCACGGTGCAGACATGCTTCAAGATGTGTATCATGTGG
 TCACCCTAAAATTCAGTTACACAGTTGCCCTAACTAGAAGACTTGCCTCCTGAACAATGGTCGCACAC
 CACAGTAAGGAATGCTCTGAAGGACTTACTGAAAGATATGAACCAGAGTTCGTTGGCCAAGGAGTGCCCC
 CTTTCACAGAGCATGATCTCCTCATTGTGAACAGCACGTAATGCAAAATGTCTCAGCAGCAAAATGTC
 AAGAATTTGGAAGGTGGTACAAACATTTCAAGAAGACAAGGATATGATGGTTGAGATGGATAGTCTGTC
 TGAACTATCCCAGCAAGGTGCCAACACGTCATTTTGGCCAGCAGCCTGTCCCAGGAAACACAGCTGAG
 CAGCCTCCATCCCCTGCGCAGCTCTCTCACGGCAGTCAGCCCTCTGTCCGGACCCCTCTTCCGAACCTGC
 ACCCTGGGCTTGTGTCAACACCGATCAGTCCTCAGCTGGTCAACCAACAGCTGGTGTAGGCTCAGTTGCT
 GAACCAGCAGTATGCAGTGAACAGACTTTAGCCAGCAGTCTTAAACCAACAGTACTTGAACCACCT
 CCCCTGTGAGTGGTCTATGAACAAGCCTTTGGAGCAGCAAGTTCCACAAACACGGAGGTCTCTTCTG
 AAATCTACCAGTGGGTGCGGGATGAACGAAACGAGCCGGAATCTCACAGCAGTATTTGCACGCGTGGC
 TTTTAACCGAACTCAGGGATTGCTTTCTGAAATCCTCCGAAAGGAAGAGGACCCCAAGACTGCATCCAG
 TCTCTGCTGGTAAACCTTCGGGCTATGAGAATTTCTTACAGTTGCCGGAAGCCGAAAGAGACCGGATAT
 ACCAGGATGAGAGGAAAGGAGCTTGAACGCAGCCTCAGCCATGGGTCTGCCCGCTGTGAGCACACC
 ACCCAGCCGCCCTCCCAGGTGAAAACAGCTACCCTTGCCACTGAGAGAAATGGGAAGCCAGAGAACAAT
 ACTATGAACATTAATGCCTCCATTTATGATGAGATTGAGCAGGAAATGAAGCGTGTCTAAAGTGTCCCAAG
 CACTGTTTGCAAAGGTTGCCGCCAAAAAGCCAGGGATGGCTGTGTGAGCTGTTGCGCTGGAAAGAAGA
 TCCTTCTCCAGAAAACAGGACCCTGTGGGAGAACCTGTGATGATCCGAAGATTCCTCAGTCTGCCTCAG
 CCTGAGCGGATGCCATCTATGAGCAGGAGAGCAATGCTGTGCATCACCATGGCGACAGACCTCCCACA
 TCATCCAGTTCCAGCAGAACAGATTACGACAGCAACAGCAGCAGCAACAGCAGCAGCAGCAGCA
 GCCACCGCCGCCACCGCCGACAGCCACAGCCACAGCCCAAGGAGGCCCCAGACTCCCCACGGCAGCCC
 ACCGTGGCCTCCTCCGCGGAGTCCGATGAGGAAAACCGGCAGAAAGACCAGGCCACGAACCAAAATTTCCG
 TGGAAGCCCTGGGGATCCTTACAGATTTTATCCAAGATGTGGGCTGTACCCAGATGAAGAGGCTATCCA
 GACTCTGTCTGCACAGCTGGACCTCCCGAAGTACACCATCATCAAGTTCTTTCAGAACCAGCGGTACTAC
 CTTAAGCACCATGGCAAGCTGAAGGACAACCTCCGGCTTGGAGGTGGATGTGGCCGAGTACAAAGATGAGG
 AGTTGCTTAAGGATTTGGAAGAGAGCGTCCAGGATAAAAAACCCCAACACCTTTTCTCAGTGAACACTAGA
 GGAAGAGCTGTGGTGAAGGGAGCACAGACGTTAATGCCGACTTGAAAGAC

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR223656 protein sequence
Red=Cloning site Green=Tags(s)

MDHLNEATQGKEHSEMSNNVSDPKGPPAKIARLEQNGSPLGRGRLGSTGGKMQGVPLKHSGLMKTNLRK
GTMLPVFCVVEHYENAI EYDCKEEHAEFVLRKDMLFNQLIEMALLSLGYSHSSAAQAKGLIQVGKWNPV
PLSYVTDAPDATVADMLQDVYHVVTLKIQLHSCPKLEDLPPEQWSHTTVRNALKDLLKDMNQSSLAKECP
LSQSMISSIVNSTYYANVSAKCOEFGRWYKHFKKTKDMMVEMDSLSEL SQQGANHVNF GQQPVPNGTAE
QPPSPAQLSHGSQPSVRTPLPNLPGLVSTPI SPQLVNQQLVMAQLLNQYAVNRLLAQQSLNQYLNHP
PPVSRSMNKPLEQVSTNTEVSSEIYQWVRDELKRAGISQAVFARVAFNRTQGLLSEILRKEEDPKTASQ
SLLVNLRAMQNFLQLPEAERDRIYQDERERSLNAASAMGPAPLLSTPPSRPPQVKATLATERNKGPENN
TMNINASIYDEIQEMKRAKVSQALFAKVAATKSQGWLCCELLRWKEDPSPENRTLWENLSMIRRFSLPQ
PERDAIYEQESNAVHHHGD RPPHI IHVPAEQIQQQQQQQQQQQQQPPPPPPQPQPQAGPRLPPRPQ
TVASSAESDEENRQKTRPRTKISVEALGILQSF IQDVGLYPDEEAIQTL SAQLDLPKYTI I KFFQNQRY
LKHHGKLDKNSGLEVDVAEYKDEELLKDLEESVQDKNANTLFSVKLEEEL SVEGSTVDNADLKD

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001163631

ORF Size: 2292 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001163631.1, NP_001157103.1</u>
RefSeq Size:	6205 bp
RefSeq ORF:	2295 bp
Locus ID:	20230
UniProt ID:	<u>Q60611</u>
Cytogenetics:	17 C
MW:	85.9 kDa
Gene Summary:	Required for the switching of fetal globin species, and beta- and gamma-globin genes regulation during erythroid differentiation. Plays a role in chromatin organization and nuclear architecture during apoptosis (By similarity). Crucial silencing factor contributing to the initiation of X inactivation mediated by Xist RNA that occurs during embryogenesis and in lymphoma. Binds to DNA at special AT-rich sequences, the consensus SATB1-binding sequence (CSBS), at nuclear matrix- or scaffold-associated regions. Thought to recognize the sugar-phosphate structure of double-stranded DNA. Transcriptional repressor controlling nuclear and viral gene expression in a phosphorylated and acetylated status-dependent manner, by binding to matrix attachment regions (MARs) of DNA and inducing a local chromatin-loop remodeling. Acts as a docking site for several chromatin remodeling enzymes and also by recruiting corepressors (HDACs) or coactivators (HATs) directly to promoters and enhancers. Modulates genes that are essential in the maturation of the immune T-cell CD8SP from thymocytes.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR223656